used as calibration standards are reanalyzed every 200 samples to verify the calibration throughout the session. If the ratio of reported to calculated concentrations for these standards drifts outside the 0.95-1.05 range during an analysis run., the cyclotron is re-tuned, system is recalibrated, and samples reanalyzed. Figure 3 shows calibration verifications and recalibrations during the analysis of March 2005 (for context) and April - June 2005 samples.

1/8 mil Mylar standards

March, April, May and June 2005 samples 1.250 1.150 Cyclotron re-tuned Cyclotron re-tuned Cvclotr ned on re-t Ratio: Found/Standard 1.050 0.950 0.850 Cyclotron re-tuned 0.750 6500 7000 7500 8000 8500 9000 9500 10000 10500 11000 11500 12000 12500 13000 13500 14000 14500 15000 Run Number

Figure 3. PESA standards for March, April, May and June 2005 samples

3.2 X-ray energy calibration

In addition to the peak counts associated with a known concentration (concentration calibration), the energy channel associated with a known fluorescence line must also be determined; this is the energy calibration. Energy calibrations were performed for the analyses of each sample month on each system to establish relationships of the form

$$energy = intercept + slope * channel$$

The following energy calibration equations (in energy units of KeV) were used for the analysis of April, May, and June 2005 samples, respectively:

- for XRF-Cu intercept= -0.05603950, slope= 0.01692728 intercept= -0.03073931, slope= 0.01665470 intercept= -0.03399334, slope= 0.01664978
- for XRF-Mo intercept= -0.07314897, slope=0.03592077 intercept=-0.07138395, slope=0.03590873 intercept=-0.07087657, slope=0.03601250

In addition, the resolution of the Si(Li) [detectors in the Mo-anode XRF systems] is frequently checked using an Fe-55 source. Results indicating changes of 5% or more in the width of the Kalpha peak for Mn are reported and further investigated. [Corrected 2/9/07.]